# <u>Standard Workflow – Non-Essential genes [ **NE1** ]</u>

Bsd first tm1 -> tm1a -> tm1/tm1a

## First Electroporation:

Diagram:			©— Bsd — 0—
Genotype:		wt/wt	tm1/wt
LOA CRIT	- OHT	2	1
LOA DEL	- OHT	2	1

tm1/wt is what you want

wt/wt is total failure of electroporation

No final loxp site with Bsd vector so no point doing + OHT assays

## Second Electroporation:

Diagram:		   •   •   •   •   •   •   •   •   •	0—1 Bsd 1—0	O Neo Neo	0 Neo Neo	O Neo Neo	o Neo Neo
Genotype:		wt/wt	tm1/wt	tm1/tm1a	wt/tm1a	tm1/tm1e	wt/tm1e
LOA CRIT	- OHT	2	1	1	2	1	1
	+ OHT	2	1	0	1	1	1
LOA DEL	- OHT	2	1	0	1	0	1
	+ OHT	2	1	0	1	0	1

tm1/tm1a is what you want

wt/tm1a is best possible control but rare

tm1/wt is second best control and very likely

wt/wt is total failure of both electroporations but still a useful control

N.B. the pattern of results here is identical to that for NE1a

## <u>Standard Workflow – Non-Essential genes [ **NE1a** ]</u>

Neo first tm1a -> tm1 -> tm1a/tm1

## First Electroporation:

Diagram:		0—————————————————————————————————————	o Neo	o Neo Neo
Genotype:		wt/wt	tm1a/wt	tm1e/wt
LOA CRIT	- OHT	2	2	2
	+ OHT	2	1	2
LOA DEL	- OHT	2	1	1
	+ OHT	2	1	1

tm1a/wt is what you want

tm1e/wt is where most of the cassette goes in but the final loxp does not

wt/wt is total failure of electroporation

## Second Electroporation:

Diagram:			o Neo	o Neo	o Neo Neo Bsd	o Neo Bsd -	0 Bsd 1 Bsd
Genotyp e:		wt/wt	tm1a/wt	tm1e/wt	tm1a/tm1	tm1e/tm1	wt/tm1
LOA CRIT	- OHT	2	2	2	1	1	1
	+ OHT	2	1	2	0	1	1
LOA DEL	- OHT	2	1	1	0	0	1
	+ OHT	2	1	1	0	0	1

tm1a/tm1 is what you want

wt/tm1 is best possible control but rare

tm1a/wt is second best control and very likely

wt/wt is total failure of both electroporations but still a useful control

N.B. the pattern of results here is identical to that for NE1

## <u>Standard Workflow – Essential genes [ E ]</u>

Neo first  $tm1a \rightarrow (Dox) \rightarrow tm1 \rightarrow tm1c/tm1$ 

## First Electroporation:

Diagram:		     		o Neo	○ Neo I	
Genotype:		wt/wt	tm1c/wt	tm1a/wt	tm1e/wt	tm1f/wt
LOA CRIT	- OHT	2	2	2	2	2
	+ OHT	2	1	1	2	2
LOA DEL	- OHT	2	1	1	1	1
	+ OHT	2	1	1	1	1
NEO		0	0	1	1	0

tm1c/wt is what you want

tm1a/wt is where the Dox does *not* excise the Neo resistance

tm1f/wt is where the Dox does excise the Neo resistance but the final loxp site was *not* inserted

tm1e/wt is where the Dox does *not* excise the Neo resistance *and* the final loxp was not inserted

wt/wt is total failure of electroporation

## <u>Standard Workflow – Essential genes [ E ]</u>

Neo first  $tm1a \rightarrow (Dox) \rightarrow tm1 \rightarrow tm1c/tm1$ 

# Second Electroporation:

Diagram:		(°	Neo Neo		o Neo		0 Bsd 1 -	o Neo Bsd -	0 Bsd -	o Neo Bsd III	0—————————————————————————————————————
Genotype:		wt/wt	tm1a/wt	tm1c/wt	tm1e/wt	tm1f/wt	wt/tm1	tm1a/tm1	tm1c/tm1	tm1e/tm1	tm1f/tm1
LOA CRIT	- OHT	2	2	2	2	2	1	1	1	1	1
	+ OHT	2	1	1	2	2	1	0	0	1	1
LOA DEL	- OHT	2	1	1	1	1	1	0	0	0	0
	+ OHT	2	1	1	1	1	1	0	0	0	0
NEO		0	1	0	1	0	0	1	0	1	0

tm1c/tm1 is what you want

tm1c/wt shows the SEP failed and is a good control

wt/tm1 is also a good control